SCIENCE CAN

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Q1: How to use a microscope stand?

A1:

1. Grooves are designed on both sides of the microscope (as shown in Figure 1), align with the slide rail on the stand (as shown in Figure 2), and slide from the top to the lowest end (as shown Figure 3).


Fig. 1


Fig. 2


Fig. 3

2. Align the LED light button on the multifunctional tool with the hole of the microscope stage and push it upward to put it inside (as shown in Figure 4).

Fig. 4

3. Rotate the focusing wheel of the stand to move the microscope upward (as shown in Figure 5), insert the specimen slide into the specimen holder and push it inward (as shown in Figure 6) until the specimen slide is stuck in the stand position (as shown in Figure 7).


Fig. 5


Fig. 6


Fig. 7

4. Rotate the focusing wheel of the stand to lower the microscope to the lowest position, press to turn on the LED light of the stand (as shown in Figure 8), and make your eyes closed to the eyepiece to start observation.


Fig. 8

5. Rotate the focusing wheel on the microscope until the image is clear (as shown Figure 9).

Fig. 9

Q2: What if I can't see anything in my field of vision?

A2:

Handheld mode:

1. Make sure the microscope LED has been turned on (as shown in Figure 10).


Fig. 10

2. Make sure that the lower end of the microscope is very much closed to the specimen slide (as shown in Figure 11).

Fig. 11

3. Confirm if the specimen is aligned with the center of the microscope objective.
4. Make observation through the eyepiece and slightly adjust the position of the specimen slide by hand until the specimen appears in the field of vision.

Stand mode:

1. Make sure the microscope LED (as shown in Figure 10) and the stand LED (as shown in Figure 12) have been turned on.

Fig. 10


Fig. 12

2. Make sure the specimen slide is stuck in the detent on the stand (as shown in Figure 13).

Fig. 13

3. Lower the microscope to the lowest point, and observe through the eyepiece to see if there is any image in the field of view. If no image can be observed, adjust the position of the observed object slightly until image can be seen.

Q3: What if the image is not at the center of the field of view?

A3:

1. Align your eye with the eyepiece to adjust the observed object slightly by hand.

2. If the image is closer to the right side of the field of view (as shown in Figure 14), move the observed object to the right (as shown in Figure 15), and the image in the field of view will move to the middle from the left.


Fig. 14

Fig. 15

3. Make adjustment in other cases using the same method. Please be noted that the actual direction of moving the observed object is opposite to the direction of movement seen in the field of vision.

Q4: What if the image is blurred?

A4:

Handheld mode:

1. Check if the lower end of the microscope is very much closed to the specimen slide (as shown in Figure 16).

Fig. 16

2. Turn the multiplier wheel to select the multiplier you need (as shown in Figure 17).

Fig. 17

3. Make your eye closed to the eyepiece, rotate the focusing wheel of the microscope, and observe the image until it becomes clear (as shown in Figure 18).

Fig. 18

Stand mode:

1. Confirm that the microscope has reached the lowest position.
2. Rotate the multiplier wheel to select the multiplier you need (as shown in Figure 19).

Fig. 19

3. Make your eye closed to the eyepiece, rotate focusing wheel of the microscope, and observe the image until it becomes clear (as shown Figure 20).


Fig. 20

Q5: How to install/remove the phone stand?

A5:

1. One side of the mobile phone stand has a pull rod that can be extended. The observation hole in the middle can slide back and forth (as shown in Figure 21).

Fig. 21

2. Prior to use, slide the observation hole to the side far away from the pull rod (as shown in Figure 22).


Fig. 22

3. Put your hand in front of the lens to confirm the main lens. The one where the hand appears on the screen is the main lens (as shown in Figure 23). Rotate the phone to the back side, to make the phone stand observation hole and the camera under the same side, and then pull the pull rod to buckle the stand on the phone (as shown in Figure 24).

Fig. 23

Fig. 24

4. Adjust the height of the mobile phone stand and the position of the observation hole slightly, to make the main lens of the camera align with the observation hole (as shown in Figure 25).

Fig. 25

5. Turn on the shooting mode of the mobile phone and confirm that the image is unobstructed (as shown in Figure 26).

Fig. 26

6. Align the threads at the inner side of the microscope eyepiece with the threads on the observation hole of the mobile phone stand, and connect the microscope and the mobile phone stand through rotation (as shown in Figure 27).

Fig. 27

7. Make observation either by holding the microscope by hand directly or placing it on the stand according to Q1 method.

8. Pull the stand pull rod to remove the mobile phone stand (as shown in Figure 28).


Fig. 28

Q6: How to use multi-function tool?

A6:

Take out/Put back the multifunctional tool body:

1. Pinch the multifunctional tool from the grooves at both sides of the bottom of the stand and take it out by pulling downward (as shown in Figure 29).
2. After using, align the LED light button on the multifunctional tool with the hole of the microscope stage (as shown in Figure 30) and put it back by pulling upward.


Fig. 29


Fig. 30

Use shovel:

1. Fasten the protruding part at the end of the shovel and rotate the shovel clockwise to take it out (as shown in Figure 31).
2. Turn the shovel clockwise to the position shown in the figure (as shown in Figure 32), and the shovel will be clamped tightly.


Fig. 31


Fig. 32

3. Use a shovel to shovel water, soil, etc. Please always be careful to prevent the microscope and toy body from contacting water to avoid short circuit problems ascribed to water.
4. After using, clean the shovel, and turn it back counterclockwise to its original position (as shown in Figure 33).


Fig. 33

Use a carving knife:

1. Fasten the protruding part of the carving knife (as shown in Figure 34) and rotate the carving knife counterclockwise to the position shown in the figure (as shown in Figure 35). At this time, the carving knife will be clamped tightly.
2. The front side of the carving knife is thin (as shown in Figure 36), which can be used for carving objects.
3. For example, thin observed objects such as leaves can be clamped in the middle of the carving knife; align with the round hole above, and press down with certain force (as shown in Figure 37) to press the observed object into small discs.
4. Clean the carving knife after using and rotate it back clockwise to its original position (as shown in Figure 38).


Fig. 34


Fig. 35


Fig. 36

Fig. 37

Fig. 38

Use tweezers:

1. Tweezers are placed at the bottom of the multifunctional tool. Buckle the protruding part of the tail end of the tweezers to take it out (as shown in Figure 39).
2. Tweezers can be used to pick up small objects. Clean the tweezers after using, and put them back in place according to the direction shown in the figure (as shown in Figure 40).


Fig. 39

Fig. 40

WARNING:


All packaging materials, such as adhesive tape, plastic sheet, binding line, labels, etc., are not part of the product and cannot be played with by children. For the safety of your children, please dispose of the above packaging materials properly. Please retain the user manual which contains important information for further reference.

This product contains a Button or Coin Cell Battery. A swallowed Button or Coin Cell Battery can cause internal chemical burns in as little as two hours and lead to death. Dispose of used batteries immediately. Keep new and used batteries away from children. If you think batteries might have been swallowed or placed inside any part of the body, seek immediate medical attention.

BATTERY INSTALLATION:

Editorial note — discrepancy in the source. Polish step 3 states “5 V AAA (LRO3/7)”, while the English text states “1.5V AAA (LRO3/AAA)” and the drawing says “AAA 1.5V”. These source statements have not been reconciled. The Polish “5 V” value is not a verified specification; confirm the correct parameter with the manufacturer or importer before replacing the battery.

1. Confirm the machine has been powered off.
2. Unscrew the screws on the battery cover of the machine (as shown in the figure).
3. Install a new 1.5V AAA (LRO3/AAA) battery with correct polarity (+, -) into the machine (as shown in the figure below).
4. Install a 3V (model CR2032) coin cell battery with the correct polarity (+, -) in the multifunction tool (as shown in the figure below).
5. Close the battery compartment cover.

PRECAUTIONS FOR BATTERY USE:

1. Alkaline battery, carbon zinc battery and other types of batteries can be used but it is forbidden to use alkaline batteries (carbon-zinc) or rechargeable batteries (nickel-cadmium) at the same time;
2. It is forbidden to recharge the non-rechargeable batteries;
3. Please recharge the rechargeable batteries under adults' supervision by only;
4. Remove the rechargeable batteries from toys before recharging;
5. It is forbidden to use different types of batteries or old and new batteries at the same time;
6. Ensure the correct polarity (+, -) of batteries;
7. Remove the depleted batteries from toys;
8. Power terminals shall not be short-circuited;
9. Please remove the battery if it will be left idled for a long time;
10. Do not throw batteries into fire;

MAINTENANCE:

1. Wipe the toy housing gently with a piece of soft, dry cloth to keep the frame clean.
2. keep the main unit away from sunlight and any other heat source.
3. Please remove the batteries if the toy is to be left unused for a long time;
4. Do not touch the main unit with any hard object or try to disassemble it.
5. Do not put the main unit in water; please keep it from moisture.


COMMON FAULTS AND TROUBLESHOOTING:

If the toy stops responding for some season,please

1. Shut down the main unit;
2. Remove the batteries and cut off power supply;
3. Keep the main unit motionless for a few min-utes, and then install the batteries again;
4. Restart the main unit and use it;
5. Please install new batteries if the main unit still doesn't work.

PRODUCT ACCESSORIES

Microscope and stand


Microscope

- Multiplier wheel
- Microscope focusing wheel
- Microscope LED light switch
- Battery cover
- Multiple display window

Microscope stand


- Stand focusing wheel
- Specimen area
- Stand LED light on/off port
- Specimen holder

Multifunctional tool and mobile phone stand

- Shovel
- Carving knife
- Tweezers
- Battery cover
- Stand LED light
- Stand LED light on/off switch
- Mobile phone stand

Other accessories

- Conventional specimen*8
- Blank specimen*4
- Label sticker*1
- Insect collection box*4
- Plant storage bag*2
- Velvet bag*1

PRODUCT ASSEMBLY

Method to assemble microscope and stand:

Align the slot on the back of the “microscope” with the slide rail at the back of the “microscope stand” and slide vertically down to the lowest position.

Method to assemble multifunction tool:

Clip the “multifunction tool” vertically upward into the bottom of the “microscope stand”. Please make sure the “raised button” faces forward and that it can snap into the hole in the “microscope stage” (as shown in the figure)

Method to assemble mobile phone stand:

1. Fix the mobile phone stand on the mobile phone, and align the center of the round hole with the center of the mobile phone camera (open the mobile phone camera to confirm that the mobile phone screen is not blocked).
2. Then screw the thread of “mobile phone stand” into the thread inside the eyepiece of “microscope”.

HOW TO USE

Handheld observation function:

Accessories required: Microscope and specimen slide

Step 1: Turn on the LED light of the “microscope”, align your eye closely with the eyepiece observation port (observer who wears glasses should take off the glasses before observation), and align the objective lens closely to the “specimen slide”.

Step 2: Adjust the position of the “specimen slide” back and forth until the specimen content appears in the field of view.

Step 3: Adjust the “multiplier wheel” at the top of the microscope (close to the eyepiece) manually until the magnification you need is reached; and then adjust the “focusing wheel” at the bottom of the microscope (close to the objective lens) until the image becomes clear.

Note: The specimen slide can be replaced by other observed objects (such as leaves, etc.). Replace it following the same operation steps.

Main observation function (specimen slide):

Accessories required: Microscope, microscope stand and specimen slide

Step 1: Assemble the “microscope” and “microscope stand”, as shown in the figure, and turn on “the stand LED light of the microscope”.

Step 2: Hold the “specimen slide” in the specimen holder until the specimen slide is clamped in place.

Step 3: Adjust the “stand focusing wheel” through rotation until the microscope tube reaches the bottom; align your eye closely to the eyepiece observation port (observer who wears glasses should take off the glasses before observation).

Step 4: Adjust the “multiplier wheel” at the top of the microscope manually (close to the eyepiece) until the magnification you need is reached; then adjust the “focusing wheel” at the bottom of the microscope (close to the objective lens) until the image becomes clear.

Main observation function (object):

Accessories required: Microscope, microscope stand and observed object (e.g. leaves)

Step 1: Assemble the “microscope” and “microscope stand” as shown in the figure, adjust the “stand focusing wheel” through rotation until the microscope tube reaches the highest position, and turn on the “the LED light of the microscope” and “stand LED light” at the same time.

Step 2: Cut the “bserved object” into dices with a diameter of about 1 cm with a carving knife and put it in the “specimen area”.

Step 3: Adjust the “stand focusing wheel” through rotation until the lower end of the “microscope tube” is close to the observed object, and align your eye closely to the observation port of the eyepiece (observer who wears glasses should take off the glasses before observation).

Step 4: Adjust the “multiplier wheel” (close to the eyepiece) at the top of the microscope through rotation until the magnification you need is reached; then adjust the “microscope focusing wheel” (close to the objective lens) until the image becomes clear.

Mobile phone observation function:

Accessories required: Microscope, mobile phone stand and specimen slide

Step 1: Assemble the “mobile phone stand” and “microscope” as shown in the figure, and turn on the microscope LED light.

Step 2: Align the microscope closely to the “specimen slide”, turn on the camera function of the mobile phone, and adjust the position of the “specimen slide” back and forth until the specimen appears in the mobile phone.

Step 3: Adjust the “multiplier wheel” (close to the eyepiece) at the top of the microscope through rotation until the magnification you need is reached; then adjust the “focusing wheel” (close to the objective lens) at the bottom of the microscope until the image becomes clear. Shoot the images for retention.

Note:

1. The specimen slide can be replaced by other observed objects (such as leaves, etc.). Replace it following the same operation steps.
2. Microscope stand can be added. The operation principle is the same. Refer to the “observation function of main body” for details about the operation steps.
3. For iPhone 12 and mobile phones above the model, camera jumps repetitively in low light environment, resulting in abnormal viewing function of the mobile phone connected with the microscope. Users are suggested to use it in a bright environment.

How to use multifunctional tool:

Accessories required: Multipurpose tool and observed object (e.g. leaves)

1. Method to use carving knife:

Rotate the knife as shown in the figure.

1. The circular area is a carving knife, and the “observed object” (note: only thin sheets such as leaves can be carved; hard objects with a thick thickness are forbidden) can be carved into a circle with a diameter of 1cm, so as to quickly make specimens.
2. The front end is a knife, which can carve the specimen into any shape required according to the specific needs.


2. Method to use shovel:

Rotate the shovel to take it out, as shown in the figure. It can be used for containing observed objects such as “water” and “soil”.

3. Method to use tweezers:

Rotate the tweezers to take them out, as shown in the figure. It can be used for clipping or tearing “plant epidermis”, “insects”, etc.

Note: Prevent the main body of multifunctional tool and microscope from contacting water during use and cleaning.

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©2025 TOPBRIGHT ANIMATION all rights reserved. MADE IN CHINA